

IOL TrueCast Guide

Prior to scanning the first time, ensure you've downloaded the latest version of the IOL TrueCast app, available for free at the App Store, and gone through the Calibration steps.

STEP 1: SCAN

Using your iPad Pro, open the IOL TrueCast app.

Click on **Scan Left** or **Scan Right**.

Position the foot so that it is illuminated in purple. Click the circular button to begin scanning, and click the circle again to end.

Spend approximately 15 seconds per foot, starting at a 45-degree angle and moving in a slight "smile" shape to capture the heel, ending at 45-degrees on the opposite side.

Review the scan. Ensure there are no holes and that the heel is well captured. If satisfied, click **Save**. Click **Rescan** if needed.


Click on **Scan Left** or **Scan Right** and repeat steps.



STEP 2: SAVE SCAN

Click the three bars in the top left hand corner, and click **Save Scan**. Choose a unique, recognizable name (patient last/first name, for example) and click **Save**.

Important Note: Naming a scan the same name as a previous scan will overwrite the previous scan, so choose a unique name.




STEP 3: OPEN THE ORDER FORM

Once you are ready to fill out your order form, click **Open Platform** to fill out using your iPad Pro.

Log in with your email and password (it will save on your devices after the first use)

Click on **New Order Form** on the top green bar.




IOL TrueCast Guide

STEP 4: FILL OUT FORM

Fill out the order form. Required fields are **Account #**, **Practitioner**, **Cast Date**, **Patient First Name**, **Patient Last Name**, **Device Type**, and **Posting Instructions**.

Once you're completed, click the top **Summary** button or **Save and Next** at the bottom. You can always navigate back to the Order Form later to make changes if necessary.



STEP 5: ATTACH SCANS

Click **Scan Mobile** on the left-hand side. Follow the prompt to open the IOL TrueCast App.

Once the app opens, click on the three bars in the top left corner, click **Open Scan List**, find the relevant scan, and click **Load Scans**. Then, click **Upload**.

Important Note: If you've the attachable mirror clip, ensure the mirror radio button is clicked "on" prior to clicking Upload.

Once the upload is confirmed successful, click **Stay on app** and swipe up from the bottom of your device to most easily navigate back to the Order Form window. Click **Refresh Scans** and green checkmarks will appear under Scan Mobile.



STEP 6: REVIEW AND SUBMIT

If you have additional images to attach, click **Add Images** from your Photo Library or Take Photo to upload up to 6 photos.

If you're submitting only one patient foot, check **Mirror Scans** to bypass the requirement for both feet.

Once satisfied, click **Submit**.

You can navigate to the Active Orders page from the browser and check on the status of your orders to ensure it has been successfully, or amend/submit in-progress orders.

